HS VFUSER

Codice Code	Volume Nominale Nominal Value	Flusso Flow Rate	Tempo Di Erogazione Nominal Infusion Time	Box
INF06120-05	60 ml	0.5 ml/h	120 h	10
INF06060-1	60 ml	1 ml/h	60 h	10
INF06030-2	60 ml	2 ml/h	30 h	10
INF06015-4	60 ml	4 ml/h	15 h	10
INF06012-5	60 ml	5 ml/h	12 h	10
INF1200-05	100 ml	0.5 ml/h	200 h	10
INF1100-1	100 ml	1 ml/h	100 h	10
INF150-2	100 ml	2 ml/h	50 h	10
INF136-28	100 ml	2.8 ml/h	36 h	10
INF125-4	100 ml	4 ml/h	25 h	10
INF120-5	100 ml	5 ml/h	20 h	10
INF117-6	100 ml	6 ml/h	17 h	10
INF113-8	100 ml	8 ml/h	13 h	10
INF110-10	100 ml	10 ml/h	10 h	10
INF3200-15	300 ml	1.5 ml/h	200 h	5
INF3150-2	300 ml	2 ml/h	150 h	5
INF360-5	300 ml	5 ml/h	60 h	5
INF338-8	300 ml	8 ml/h	38 h	5
INF330-10	300 ml	10 ml/h	30 h	5
DISPOSITIVI DEDICATI ALL'INFUSIONE DI BREVE DURATA (ALTO FLUSSO) MODELS FOR SHORT THERAPIES (HIGH FLOW)				
INF 051-50	50 ml	50 ml/h	1 h	10
INF 12-50	100 ml	50 ml/h	2 h	10
INF 11-100	100 ml	100 ml/h	1 h	10
INF 255-50	250 ml	50 ml/h	5 h	5
INF 252.5-100	250 ml	100 ml/h	2.5 h	5
INF 251-250	250 ml	250 ml/h	1 h	5
DISPOSITIVI DEDICATI ALL'INFUSIONE DI FLUOROURACILE MODELS FOR FLUOROURACIL INFUSION				
INF1200-05 FU	100 ml	0.5 ml/h	200 h	10
INF150-2FU	100 ml	2 ml/h	50 h	10

HS HOSPITAL SERVICE S.P.A
Sede legale/Registered office:
Via Zosimo 13 - 00178 Roma, ITALY

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Pompe Elastomeriche Infusionali
per terapie ospedaliere o domiciliari

Elastomeric Infusional Pumps
for hospital or home treatments

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ELASTOMERIC PUMP FOR CONTINUOUS INFUSION WITH FIXED FLOW RATE

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HS VFUSER is a device for continuous intra-venous, intra-arterial, subcutaneous, epidural, subarachnoid or peri-nervous infusion of drugs (analgesics, anaesthetics, chemotherapeuticals, antibiotics and so on).

Device features:

- silicone balloon reservoir
- UV protected transparent, rigid and graduated cover to check the balloon
- UV protected infusion line with star-shaped section, anti-kinking and twisting resistant, already connected to the balloon
- one-way injection valve with luer lock connector and cap
- inner 5um anti-particle filter
- latex and phthalates free
- color-coded label and indication of flow rate on the device
- variable filling: the device may be filled partially and not necessarily at the nominal volume
- accuracy $\pm 10\%$
- the maximum tolerance over the nominal volume is +10%
- traceability: the manufacturing lot number is etched on the device

Accessories: transport bag, patient's data label, traceability label, three way stopcock, luer lock cap (for FU models a safety connector is included).

If the drug Fluorouracil must be administered, we recommend to use dedicated HS VFUSER (codes ending with FU).

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ELASTOMERIC PUMP FOR CONTINUOUS INFUSION WITH VARIABLE FLOW RATE

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MULTIFLOW is a device for continuous intra-venous, intra-arterial, subcutaneous, epidural, subarachnoid or peri-nervous infusion of drugs (analgesics, anaesthetics, chemotherapeuticals, antibiotics and so on).

Device features:

- silicone balloon reservoir
- UV protected transparent, rigid and graduated cover to check the balloon
- UV protected infusion line with star-shaped section, anti-kinking and twisting resistant, already connected to the balloon
- one-way injection valve with luer lock connector and cap
- inner 5um anti-particle filter
- selector to choose the flow-rate between the 11 available flow rates, with safety key to lock the chosen flow rate and to avoid accidental handlings or dislocations
- latex and phthalates free
- color-coded label and indication of flow rate on the device
- variable filling: the device may be filled partially and not necessarily at the nominal volume
- accuracy $\pm 10\%$
- the maximum tolerance over the nominal volume is +10%
- traceability: the manufacturing lot number is etched on the device

Accessories: transport bag, patient's data label, traceability label, three way stopcock, luer lock cap.

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ELASTOMERIC PUMP FOR CONTINUOUS INFUSION WITH FIXED FLOW RATE WITH PCA (PATIENT-CONTROLLED ANALGESIA) MODULE

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HS VFUSER PCA is a device for continuous infusion of drugs with PCA module for the additional administration of the drug at the patient request.

Device features:

- silicone balloon reservoir
- UV protected transparent, rigid and graduated cover to check the balloon
- UV protected infusion line with star-shaped section, anti-kinking and twisting resistant, already connected to the balloon
- one-way injection valve with luer lock connector and cap
- inner 5um anti-particle filter
- PCA module assembled to the infusion line, with waiting time between subsequent deliveries of boluses of 30 minutes
- latex and phthalates free
- color-coded label and indication of flow rate on the device
- variable filling: the device may be filled partially and not necessarily at the nominal volume
- accuracy $\pm 10\%$
- the maximum tolerance over the nominal volume is +10%
- traceability: the manufacturing lot number is etched on the device

Accessories: transport bag, patient's data label, traceability label, three way stopcock, luer lock cap.

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ELASTOMERIC PUMP FOR CONTINUOUS INFUSION WITH VARIABLE FLOW RATE WITH PCA (PATIENT-CONTROLLED ANALGESIA) MODULE

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MULTIFLOW PCA is a device for continuous infusion of drugs with PCA module for the additional administration of the drug at the patient request.

Device features:

- silicone balloon reservoir
- UV protected transparent, rigid and graduated cover to check the balloon
- UV protected infusion line with star-shaped section, anti-kinking and twisting resistant, already connected to the balloon
- one-way injection valve with luer lock connector and cap
- inner 5um anti-particle filter
- selector to choose the flow-rate between the 11 available flow rates, with safety key to lock the chosen flow rate and to avoid accidental handlings or dislocations
- PCA module assembled to the infusion line, with waiting time between subsequent deliveries of boluses of 30 minutes
- latex and phthalates free
- color-coded label and indication of flow rate on the device
- variable filling: the device may be filled partially and not necessarily at the nominal volume
- accuracy $\pm 10\%$
- the maximum tolerance over the nominal volume is +10%
- traceability: the manufacturing lot number is etched on the device

Accessories: transport bag, patient's data label, traceability label, three way stopcock, luer lock cap.