
AMICA-GEN

PROGRAMMABLE GENERATOR
FOR MICROWAVE AND/OR RADIOFREQUENCY
ENERGY DELIVERY AND CONTROL

1 SYSTEM 2 TECHNOLOGIES MWA & RFA

*FREEDOM OF CHOICE
HIGH PERFORMANCE
OPTIMIZED CONTROL
ENHANCED SAFETY
EXTENDED INTERACTIVITY
DEDICATED SOFTWARE*

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PRODUCT DESCRIPTION

AMICA-GEN is a programmable microwave and/or radiofrequency generator for thermal ablation treatments.

- *Easy and user-friendly graphic interface with touch screen display to set-up the procedure parameters*
- *Display or footswitch-controlled energy activation/shutdown*
- *Peristaltic pump built into the generator with automatic control*
- *Automatic identification of the applicator and selection of the corresponding default working parameters*
- *Separate ports for the connection and monitoring of the grounding pads (for models comprising radiofrequency module)*
- *Continuous monitoring of the relevant parameters during the ablation: time, delivered power, temperature, reflected power (only in microwave treatments), impedance (only in radiofrequency treatments), energy deposited in the patient*
- *Maximum portability and compactness. Dimensions: 47 cm x 35 cm x 13 cm*

- **MW energy deployment modalities:**
 - **MAN** mode: continuous delivery of the set power (up to 80W net on probe)
 - **PULSE** mode: intermittent delivery of the set power (up to 180W net on probe -205W for -S1.5 probes, with short cables -)

- **RF energy deployment modalities:**
 - **MAN** mode: continuous delivery of the set power and visual/acoustic notification of impedance hikes
 - **AUTO** mode: automatic power adjustment (up to the set power), based on impedance changes
 - **TEMP** mode: automatic adjustment of the output power for reaching and maintaining a manually selected tissue temperature target

- *Automatic RF and MW track ablation*

- *Sequential or simultaneous multi-probe operation on 2 or 3 MW antennas or RF electrodes through the **AMICA-SWITCHBOX** accessory, also usable as thermocouple reader (with 3 measurement ports) for perilesional thermal monitoring.*

* **AMICA-SWITCHBOX** may be used only in conjunction with Amica-GEN units featuring firmware v1.8.0 or later

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* Available accessories for moving and standing of AMICA-GEN not part of HS AMICA system.