THYROID ABLATION

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INTRODUCTION

Thyroid nodules are an extremely common occurrence. Until recently, international guidelines suggested either no treatment or surgery for patients with benign solid thyroid nodules, depending on their size, cytology, and symptoms.

Lately, several non-surgical, minimally invasive approaches for treating thyroid nodules have been developed and, nowadays, they provide a valid alternative to surgery for treating symptomatic benign thyroid nodules. These techniques include percutaneous ethanol injection therapy (PEIT), laser ablation therapy (LAT), and radiofrequency ablation (RFA) or microwave ablation (MWA).

Among these techniques, RFA seems to provide an optimal trade-off between several requirements, such as treatment efficacy, safety, adaptability to different clinical needs (e.g. recurrent cysts, mixed or solid benign thyroid nodules, hyperthyroidism), cosmetic compliance, intra-operative and post-operative pain management and cost-effectiveness.

However, the most recent technological advancements in MWA increased its safety to the same level of RFA treatments, but with a remarkably higher heating velocity, that proves beneficial when dealing with large, or hypervascularized, or cystic nodules.

INDICATIONS

- *Presence of symptoms (neck pain, dysphagia, foreign body sensation, discomfort and cough)*
- *Cosmetic concerns*
- *AFTN (Autonomously Functioning Thyroid Nodules)*
- *Recurrent cancer (in pts. at high surgical risk or pts. who refuse repeated surgery)*
- ***Therapeutic goal will be VOLUME Reduction***

MOVING-SHOT TECHNIQUE

The most used treatment technique is the so-called “moving-shot” technique, which consists in performing multiple overlapping ablations in contiguous positions within the target nodule (Fig. 1): this approach is safe, effective and easily conformable to any nodule size and shape. HS AMICA features dedicated RF and MW energy delivery modes for facilitating the deployment of moving shot treatments..


Fig.1: Sketch of the moving-shot technique for RFA treatment of benign thyroid nodules

In both cases, the HS AMICA system detects and records the status of target tissues at the beginning of each moving shot and keeps track of any change occurring during the ablation: the progressive dehydration induced by the thermal treatment causes the impedance measured between the RF electrode tip and the grounding pads or, respectively, the reflection coefficient measured by the MW antenna to rise.

The increased value of these quantities beyond a given threshold suggests that the coagulation of the tissues closest to the RF or MW probe distal tip has been achieved. If this is the case, HS AMICA notifies the user through an acoustic signal and automatically reduces the power output. Normal energy delivery at the power level set by the user is automatically reprised -i.e. a new moving shot starts- upon moving of the probe tip to tissues not yet coagulated, as indicated by the return of RF impedance or MW reflection coefficient measurements to their baseline values.

On the RF/MW generator the moving shot counter and the deposited energy counter are constantly updated.

MOVING SHOT TREATMENTS WITH RF AMICA-PROBE: THE RIMA APPROACH

HS AMICA provides for a dedicated RF energy delivery algorithm for the semi-automated handling of moving shot treatments, based on the monitoring of the variations in tissue impedance (**RIMA: RF IMPEDANCE-BASED MOVING SHOT ALGORITHM**).

The RF Amica-Probe models for which RIMA mode is available are those listed in the table below.

Code	Description
RFH18070E05V1*	18G x 70mm - exposed tip 5mm - 90° handle
RFH18070E07V1*	18G x 70mm - exposed tip 7mm - 90° handle
RFH18070E10V1*	18G x 70mm - exposed tip 10mm - 90° handle
RFH18070E15V1*	18G x 70mm - exposed tip 15mm - 90° handle
RFH18100E05V1*	18G x 100mm - exposed tip 5mm - 90° handle
RFH18100E07V1*	18G x 100mm - exposed tip 7mm - 90° handle
RFH18100E10V1*	18G x 100mm - exposed tip 10mm - 90° handle
RFH18100E15V1*	18G x 100mm - exposed tip 15mm - 90° handle

* For the version with straight (pencil-style) handle, replace V1 with V2

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HS AMICA provides for a dedicated MWenergy delivery algorithm for the semi-automated handling of moving shot treatments, based on the monitoring of the variations in the reflection coefficient measured by the MW antenna (**MONTREAL: MICROWAVE MOVING SHOT REFLECTION-BASED ALGORITHM**).

The Amica-Probe models for which MONTREAL mode is available are those listed in the table below.

Code	Description
APK17070T19V5	17G x 7cm - 90° handle
APK17100T19V5	17G x 10cm - 90° handle
APK17150T19V5	17G x 15cm - 90° handle
APK17070T19V6	17G x 7cm - straight handle
APK17100T19V6	17G x 10cm - straight handle
APK17150T19V6	17G x 15cm - straight handle
APK18070T19V5	18G x 7cm - 90° handle
APK18100T19V5	18G x 10cm - 90° handle
APK18150T19V5	18G x 15cm - 90° handle
APK18070T19V6	18G x 7cm - straight handle
APK18100T19V6	18G x 10cm - straight handle
APK18150T19V6	18G x 15cm - straight handle