

ALM SET

REPERE PER LA LOCALIZZAZIONE DELLE LESIONI DELLA MAMMELLA BREAST LOCALIZATION NEEDLE

ALM SET è un ago con repere non riposizionabile per la localizzazione pre-operatoria delle lesioni non palpabili della mammella.

Caratteristiche del dispositivo:


- ago centimetrato, con punta lancet, fermo di profondità e connettore luer lock con invito ad imbuto
- repere metallico non riposizionabile con punta a V o Z e marker di riferimento per l'individuazione della posizione della punta del repere

ALM SET is a needle with a non-repositionable wire intended for pre-surgical localization of non-palpable breast lesions.

Device Features:

- needle with lancet tip, depth marks, depth stop and throated luer lock connector for easy wire loading
- non-repositionable metallic wire with V or Z tip and reference marks to identify anchor position

CODE	DIAMETER	LENGTH	BOX
ALM SET Z TYPE – punta a Z/Z type tip			
ALMZ 2005	20G	50mm	20
ALMZ 20075	20G	75mm	20
ALMZ 2010	20G	100mm	20
ALMZ 2015	20G	150mm	20
ALM SET V TYPE – punta a V/V type tip			
ALM 2005	20G	50mm	20
ALM 20075	20G	75mm	20
ALM 2010	20G	100mm	20
ALM 2015	20G	150mm	20
ALM 2105	21G	50mm	20
ALM 21075	21G	75mm	20
ALM 2110	21G	100mm	20



PUNTA DEL REPERE V o Z TYPE
V o Z TYPE WIRE TIP

