

# BMN J TYPE

## SET PER BIOPSIA OSTEOMIDOLLARE BONE MARROW BIOPSY SET

**BMN J TYPE** è un ago per biopsia osteomidollare tipo Jamshidi.

### Caratteristiche del dispositivo:

- ago tipo Jamshidi rastremato, centimetrato, con punta a bocca di pesce con impugnatura ergonomica con connettore luer lock che consente di distribuire la forza su tutto il palmo della mano e non localmente
- mandrino estraibile, perfettamente accoppiato all'ago, con punta piramidale a tre facce
- filo di estrazione a punta smussa, centimetrato, indicato per il recupero del campione
- tappo di protezione ("safety") indicato per ricoprire la punta dell'ago ed evitare pericoli per l'operatore
- tappo luer lock indicato per chiudere l'ago evitando la fuoriuscita di sangue

**BMN J TYPE** is a Jamshidi type bone marrow biopsy needle.

### Device Features:

- Jamshidi needle with tapered "fish mouth" tip, depth marks and ergonomic handle with luer lock connector (it ensures safe grip and manoeuvrability and force distribution on the palm of the hand)
- removable stylet with pyramidal tip
- extraction wire with blunt tip intended for sample retrieval
- safety cap intended for needle tip covering to avoid damages for personnel
- luer lock cap

| CODE   | DIAMETER | LENGTH | BOX |
|--------|----------|--------|-----|
| J 0810 | 8G       | 100mm  | 10  |
| J 0815 | 8G       | 150mm  | 10  |
| J 1110 | 11G      | 100mm  | 10  |

| CODE   | DIAMETER | LENGTH | BOX |
|--------|----------|--------|-----|
| J 1115 | 11G      | 150mm  | 10  |
| J 1306 | 13G      | 60mm   | 10  |
| J 1310 | 13G      | 100mm  | 10  |

FUTUREMADEPRESENT

