

DISPOSITIVO SEMIAUTOMATICO PER BIOPSIA ISTOLOGICA CON TECNICA MENGHINI MODIFICATA SEMI-AUTOMATIC DEVICE FOR MODIFIED MENGHINI-TYPE BIOPSY

BIOMOL è un ago per biopsia dei tessuti molli semiautomatico che automatizza la procedura, lasciando allo stesso tempo all'operatore la sensibilità del prelievo.

Il BIOMOL permette di eseguire la procedura biottica con una sola mano, rendendola più veloce rispetto agli aghi manuali.

Caratteristiche del dispositivo:

- cannula sganciabile con punta Menghini, centimetrata
- mandrino a punta piramidale solidale al pistone della siringa
- depressione automatica
- pulsante di azionamento laterale
- marker ecogeno interno

BIOMOL is a soft tissues biopsy needle. It automates the biopsy procedure and it leaves full control of the sampling to the physician at the same time.


The BIOMOL allows to perform the biopsy with one hand only, making it faster versus manual devices.

Needle features:

- removable cannula with Menghini tip and depth marks
- plunger integrated stylet with pyramidal tip
- automatic depression block
- lateral firing button
- inner echogenic marker

CODE	DIAMETER	LENGTH	BOX
BM1607	16G	70mm	10
BM1609	16G	90mm	10
BM1612	16G	120mm	10
BM1615	16G	150mm	10
BM1620	16G	200mm	10
BM1709	17G	90mm	10
BM1715	17G	150mm	10
BM1720	17G	200mm	10
BM1807	18G	70mm	10
BM1810	18G	100mm	10
BM1815	18G	150mm	10
BM1820	18G	200mm	10
BM1910	19G	100mm	10

CODE	DIAMETER	LENGTH	BOX
BM1915	19G	150mm	10
BM1920	19G	200mm	10
BM2010	20G	100mm	10
BM2015	20G	150mm	10
BM2020	20G	200mm	10
BM2110	21G	100mm	10
BM2115	21G	150mm	10
BM2120	21G	200mm	10



PULSANTE LATERALE PER
DEPRESSIONE AUTOMATICA
LATERAL PUSH-BUTTON FOR
AUTOMATIC DEPRESSION



MENGHINI Type



HS[®]
FUTUREMADEPRESENT