

ECOJEKT® DIVA

AGO PER BIOPSIA CITOLOGICA ED ALTRE APPLICAZIONI
NEEDLE FOR FNAC AND OTHER APPLICATIONS

ECOJEKT DIVA è un ago con punta tipo Chiba a diametro variabile. Esso è caratterizzato da una variazione di 2G tra il diametro della punta e quello del corpo della cannula. Ciò consente di pungere il bersaglio con un diametro piccolo ma con la garanzia di stabilità e precisione.

Caratteristiche del dispositivo:

- ago con punta a doppia affilatura, centimetrato, con fermo di profondità
- mandrino rimovibile con codice colore
- marker ecogeno interno
- applicazioni: biopsia citologica, agoaspirazione, amniocentesi, villocentesi, prelievo di ovociti, prelievo di sangue fetale, alcolizzazione, puntura esplorativa, iniezione di mezzi di contrasto o altre sostanze

ECOJEKT DIVA is a needle with Chiba type tip and variable diameter. It features a 2G difference between the tip diameter and the cannula diameter, allowing to reach the target with a little diameter and guaranteeing stability and precision.

Device Features:

- needle with double sharpening tip, depth marks and depth stop
- removable colour coded stylet
- inner echogenic marker
- intended uses: cytological biopsy, fine-needle aspiration, amniocentesis, chorionic villus sampling, oocyte pick-up, fetal blood sampling, alcoholization, exploratory puncture, contrast medium and other fluids injection

CODE	DIAMETER	LENGTH	BOX
DIVA2010	18/20G	100mm	20
DIVA2015	18/20G	150mm	20
DIVA2020	18/20G	200mm	20
DIVA2107	19/21G	70mm	20
DIVA2110	19/21G	100mm	20
DIVA2115	19/21G	150mm	20
DIVA2120	19/21G	200mm	20
DIVA2210	20/22G	100mm	20
DIVA2215	20/22G	150mm	20

CODE	DIAMETER	LENGTH	BOX
DIVA2220	20/22G	200mm	20
DIVA2307	21/23G	70mm	20
DIVA2310	21/23G	100mm	20
DIVA2315	21/23G	150mm	20
DIVA2320	21/23G	200mm	20
DIVA2407	22/24G	70mm	20
DIVA2410	22/24G	100mm	20
DIVA2415	22/24G	150mm	20
DIVA2608	24/26G	80mm	20



DETTAGLIO DELLA PUNTA
TIP DETAIL

