

HS BILIARY

CATETERE PER DRENAGGIO BILIARE ESTERNO/INTERNO NEPHROSTOMY DRAINAGE CATHETER

HS BILIARY è un catetere per drenaggio biliare esterno/interno.

Caratteristiche del dispositivo:

- catetere in poliuretano, flessibile, radiopaco, con rivestimento idrofilico, resistente al kinking e all'ostruzione
- punta rinforzata per evitare l'incamiciamento durante l'inserimento
- punta pigtail (2cm) rastremata con ampi fori
- fori addizionali sul corpo del catetere
- disponibili modelli con sistema di ancoraggio del pigtail
- connettore luer lock
- include raddrizzatore per la punta, cannula di irrigidimento metallica a punta smussa, cannula flessibile e tappo/valvola
- eccellente capacità di drenaggio

HS BILIARY is an external/internal biliary drainage catheter.

Catheter features:

- polyurethane catheter, flexible, radiopaque, with hydrophilic coating, kink and encrustation resistant
- reinforced tip to help prevent the catheter from buckling during insertion
- 2cm pigtail tip with large holes
- additional holes on the catheter
- "locking system" models available
- luer lock connector
- includes split straightener, metal cannula, flexible cannula and septum cap
- excellent flow rates

CODE	DIAMETER	WORKING LENGTH	BOX
Non-locking catheters			
HSBP0850	Fr. 8	490mm	5
HSBP1050	Fr. 10	490mm	5
HSBP1250	Fr. 12	490mm	5
HSBP1450	Fr. 14	490mm	5
Locking catheters			
HSBPL0850	Fr. 8	490mm	5
HSBPL1040	Fr. 10	390mm	5
HSBPL1050	Fr. 10	490mm	5
HSBPL1250	Fr. 12	490mm	5
HSBPL1450	Fr. 14	490mm	5

*I cateteri Fr.8-14 accettano guide 0.038". Altre lunghezze disponibili su richiesta.
Fr.8-14 catheters accept 0.038" guidewire. Other lengths available upon request.*

FUTUREMADEPRESENT

