

HS GENERAL PURPOSE MINIPIGTAIL

CATETERE PER DRENAGGIO PERCUTANEO PERCUTANEOUS DRAINAGE CATHETER

HS GENERAL PURPOSE MINIPIGTAIL è un catetere per uso generale con tecnica diretta.

Caratteristiche del dispositivo:

- catetere in poliuretano, flessibile, radiopaco, con rivestimento idrofilico, resistente al kinking e all'ostruzione
- punta rinforzata per evitare l'incamiciamento durante l'inserimento
- punta pigtail (1cm) rastremata con ampi fori
- sistema di ancoraggio del pigtail
- connettore luer lock
- include raddrizzatore per la punta, mandrino a punta piramidale e cannula di irrigidimento metallica a punta smussa
- eccellente capacità di drenaggio

HS GENERAL PURPOSE MINIPIGTAIL is a general purpose catheter intended to be used with one-shot technique.

Device Features:

- polyurethane catheter, flexible, radiopaque, with hydrophilic coating, kink and encrustation resistant
- reinforced tip to help prevent the catheter from buckling during insertion
- 1cm pigtail tip with large holes
- locking system for pigtail tip
- luer lock connector
- includes split straightener, trocar and metal cannula
- excellent flow rates

CODE	DIAMETER	WORKING LENGTH	BOX
HSGPML0620	Fr. 6	190mm	5
HSGPML0625	Fr. 6	240mm	5
HSGPML0830	Fr. 8	290mm	5
HSGPML1030	Fr. 10	290mm	5

*I cateteri Fr.6 accettano guide 0.035". I cateteri Fr.8-10 accettano guide 0.038".
Fr.6 catheters accept 0.035" guidewire. Fr.8-10 catheters accept 0.038" guidewire.*

FUTUREMADEPRESENT

