

HS GENERAL PURPOSE

CATETERE PER DRENAGGIO PERCUTANEO PERCUTANEOUS DRAINAGE CATHETER

HS GENERAL PURPOSE è un catetere per uso generale con tecnica diretta o Seldinger.

Caratteristiche del dispositivo:

- catetere in poliuretano, flessibile, radiopaco, con rivestimento idrofilico, resistente al kinking e all'ostruzione
- punta rinforzata per evitare l'incamiciamento durante l'inserimento
- punta pigtail (2cm) rastremata con ampi fori
- disponibili modelli con sistema di ancoraggio del pigtail
- connettore luer lock
- include raddrizzatore per la punta, mandrino a punta piramidale e cannula di irrigidimento metallica a punta smussa
- eccellente capacità di drenaggio

HS GENERAL PURPOSE is a general purpose catheter intended to be used with one-shot or Seldinger technique.

Device Features:

- polyurethane catheter, flexible, radiopaque, with hydrophilic coating, kink and encrustation resistant
- reinforced tip to help prevent the catheter from buckling during insertion
- 2cm pigtail tip with large holes
- "locking system" models available
- luer lock connector
- includes split straightener, trocar and metal cannula
- excellent flow rates

CODE	DIAMETER	WORKING LENGTH	BOX
Non-locking catheters			
HSGP0530	Fr. 5	290mm	5
HSGP0625	Fr. 6	240mm	5
HSGP0630	Fr. 6	290mm	5
HSGP0725	Fr. 7	240mm	5
HSGP0730	Fr. 7	290mm	5
HSGP0740	Fr. 7	390mm	5
HSGP0825	Fr. 8	240mm	5
HSGP0830	Fr. 8	290mm	5
HSGP0840	Fr. 8	390mm	5
HSGP1030	Fr. 10	290mm	5
HSGP1040	Fr. 10	390mm	5
HSGP1230	Fr. 12	290mm	5
HSGP1240	Fr. 12	390mm	5
HSGP1430	Fr. 14	290mm	5
HSGP1440	Fr. 14	390mm	5
HSGP1630	Fr. 16	290mm	5
HSGP1640	Fr. 16	390mm	5

CODE	DIAMETER	WORKING LENGTH	BOX
Locking catheters			
HSGPL0625	Fr. 6	240mm	5
HSGPL0630	Fr. 6	290mm	5
HSGPL0725	Fr. 7	240mm	5
HSGPL0730	Fr. 7	290mm	5
HSGPL0740	Fr. 7	390mm	5
HSGPL0825	Fr. 8	240mm	5
HSGPL0830	Fr. 8	290mm	5
HSGPL0840	Fr. 8	390mm	5
HSGPL1030	Fr. 10	290mm	5
HSGPL1040	Fr. 10	390mm	5
HSGPL1230	Fr. 12	290mm	5
HSGPL1240	Fr. 12	390mm	5
HSGPL1430	Fr. 14	290mm	5
HSGPL1440	Fr. 14	390mm	5
HSGPL1630	Fr. 16	290mm	5
HSGPL1640	Fr. 16	390mm	5

I cateteri Fr.6-7 accettano guide 0.035". I cateteri Fr.8-16 accettano guide 0.038". Altre lunghezze disponibili su richiesta.
Fr.6-7 catheters accept 0.035" guidewire. Fr.8-16 catheters accept 0.038" guidewire. Other lengths available upon request.

FUTUREMADEPRESENT

