

HS NEPHRO CATHETER KIT

KIT CON CATETERE PER DRENAGGIO NEFROSTOMICO
NEPHROSTOMY DRAINAGE CATHETER KIT

HS NEPHRO CATHETER KIT è un kit per drenaggio nefrostomico.

Caratteristiche del dispositivo:

- catetere in poliuretano, flessibile, radiopaco, con rivestimento idrofilico, resistente al kinking e all'ostruzione
- punta rinforzata per evitare l'incamiciamento durante l'inserimento
- punta pigtail (2cm) rastremata con ampi fori
- disponibili modelli con sistema di ancoraggio del pigtail
- connettore luer lock
- include raddrizzatore per la punta, cannula di irrigidimento metallica a punta smussa e cannula flessibile
- eccellente capacità di drenaggio

HS NEPHRO CATHETER KIT include:

ago introduttore 18G, guida flessibile o tipo Lunderquist (teflonata, idrofilica, con tratto terminale floppy e punta a "J" 3mm), dilatatori sequenziali con punta rastremata, raccordo di estensione con connettori luer lock, raccordo per sacca con connettori luer lock e ad imbuto, ago 22G per iniezione mezzo di contrasto, dispositivo di stabilizzazione, valvola per il lavaggio del catetere, telo per campo sterile, salvietta bianca, garze (5), rubinetto a tre vie, bisturi, applicatore con spugna.

Nella confezione è inclusa una bustina di gel sterile da 20g Aquasonic 100.

HS NEPHRO CATHETER KIT is a nephrostomy drainage catheter.

Device Features:

- *polyurethane catheter, flexible, radiopaque, with hydrophilic coating, kink and encrustation resistant*
- *reinforced tip to help prevent the catheter from buckling during insertion*
- *2cm pigtail tip with large holes*
- *"locking system" models available*
- *luer lock connector*
- *includes split straightener, metal cannula and flexible cannula*
- *excellent flow rates*

HS NEPHRO CATHETER KIT includes:

18G introducer needle, flexible or Lunderquist type guidewire (PTFE-coated, hydrophilic, with floppy and 3mm "J" tip), dilators, luer lock extension tube, funnel-luer lock adapter, 22G needle for contrast medium injection, catheter stabilizer, septum cap, drape, white drape, gauzes (5), three-way stopcock, scalpel, sponge applicator.

Into the box a 20g Aquasonic 100 gel sterile packet is included.

CODE	DIAMETER	WORKING LENGTH	BOX
<i>Non-locking catheters and flexible guidewire</i>			
HSNP0740KF	Fr. 7	390mm	1
HSNP0750KF	Fr. 7	490mm	1
HSNP0840KF	Fr. 8	390mm	1
HSNP0850KF	Fr. 8	490mm	1
HSNP1040KF	Fr. 10	390mm	1
HSNP1050KF	Fr. 10	490mm	1
HSNP1240KF	Fr. 12	390mm	1
HSNP1250KF	Fr. 12	490mm	1
HSNP1440KF	Fr. 14	390mm	1
HSNP1450KF	Fr. 14	490mm	1

CODE	DIAMETER	WORKING LENGTH	BOX
<i>Non-locking catheters and Lunderquist guidewire</i>			
HSNP0740KL	Fr. 7	390mm	1
HSNP0750KL	Fr. 7	490mm	1
HSNP0840KL	Fr. 8	390mm	1
HSNP0850KL	Fr. 8	490mm	1
HSNP1040KL	Fr. 10	390mm	1
HSNP1050KL	Fr. 10	490mm	1
HSNP1240KL	Fr. 12	390mm	1
HSNP1250KL	Fr. 12	490mm	1
HSNP1440KL	Fr. 14	390mm	1
HSNP1450KL	Fr. 14	490mm	1

CODE	DIAMETER	WORKING LENGTH	BOX
<i>Locking catheters and flexible guidewire</i>			
HSNPL0740KF	Fr. 7	390mm	1
HSNPL0750KF	Fr. 7	490mm	1
HSNPL0840KF	Fr. 8	390mm	1
HSNPL0850KF	Fr. 8	490mm	1
HSNPL1040KF	Fr. 10	390mm	1
HSNPL1050KF	Fr. 10	490mm	1
HSNPL1240KF	Fr. 12	390mm	1
HSNPL1250KF	Fr. 12	490mm	1
HSNPL1440KF	Fr. 14	390mm	1
HSNPL1450KF	Fr. 14	490mm	1

CODE	DIAMETER	WORKING LENGTH	BOX
<i>Locking catheters and Lunderquist guidewire</i>			
HSNPL0740KL	Fr. 7	390mm	1
HSNPL0750KL	Fr. 7	490mm	1
HSNPL0840KL	Fr. 8	390mm	1
HSNPL0850KL	Fr. 8	490mm	1
HSNPL1040KL	Fr. 10	390mm	1
HSNPL1050KL	Fr. 10	490mm	1
HSNPL1240KL	Fr. 12	390mm	1
HSNPL1250KL	Fr. 12	490mm	1
HSNPL1440KL	Fr. 14	390mm	1
HSNPL1450KL	Fr. 14	490mm	1