

HS NEPHRO CATHETER

CATETERE PER DRENAGGIO NEFROSTOMICO NEPHROSTOMY DRAINAGE CATHETER

HS NEPHRO CATHETER è un catetere per drenaggio nefrostomico.

HS NEPHRO CATHETER is a nephrostomy drainage catheter.

Caratteristiche del dispositivo:

- catetere in poliuretano, flessibile, radiopaco, con rivestimento idrofilico, resistente al kinking e all'ostruzione
- punta rinforzata per evitare l'incamiciamento durante l'inserimento
- punta pigtail (2cm) rastremata con ampi fori
- disponibili modelli con sistema di ancoraggio del pigtail
- connettore luer lock
- include raddrizzatore per la punta, cannula di irrigidimento metallica a punta smussa e cannula flessibile
- eccellente capacità di drenaggio

Device Features:

- polyurethane catheter, flexible, radiopaque, with hydrophilic coating, kink and encrustation resistant
- reinforced tip to help prevent the catheter from buckling during insertion
- 2cm pigtail tip with large holes
- "locking system" models available
- luer lock connector
- includes split straightener, metal cannula and flexible cannula
- excellent flow rates

CODE	DIAMETER	WORKING LENGTH	BOX
Non-locking catheters			
HSNP0740	Fr. 7	390mm	5
HSNP0750	Fr. 7	490mm	5
HSNP0840	Fr. 8	390mm	5
HSNP0850	Fr. 8	490mm	5
HSNP1040	Fr. 10	390mm	5
HSNP1050	Fr. 10	490mm	5
HSNP1240	Fr. 12	390mm	5
HSNP1250	Fr. 12	490mm	5
HSNP1440	Fr. 14	390mm	5
HSNP1450	Fr. 14	490mm	5

CODE	DIAMETER	WORKING LENGTH	BOX
Locking catheters			
HSNPL0630	Fr. 6	290mm	5
HSNPL0730	Fr. 7	290mm	5
HSNPL0740	Fr. 7	390mm	5
HSNPL0750	Fr. 7	490mm	5
HSNPL0830	Fr. 8	290mm	5
HSNPL0840	Fr. 8	390mm	5
HSNPL0850	Fr. 8	490mm	5
HSNPL1030	Fr. 10	290mm	5
HSNPL1040	Fr. 10	390mm	5
HSNPL1050	Fr. 10	490mm	5
HSNPL1230	Fr. 12	290mm	5
HSNPL1240	Fr. 12	390mm	5
HSNPL1250	Fr. 12	490mm	5
HSNPL1440	Fr. 14	390mm	5
HSNPL1450	Fr. 14	490mm	5

I cateteri Fr.7 accettano guide 0.035". I cateteri Fr.8-14 accettano guide 0.038". Altre lunghezze disponibili su richiesta.
Fr.7 catheters accept 0.035" guidewire. Fr.8-14 catheters accept 0.038" guidewire. Other lengths available upon request.

FUTUREMADEPRESENT

