

HS NEPHRO MALECOT

CATETERE PER DRENAGGIO NEFROSTOMICO NEPHROSTOMY DRAINAGE CATHETER

HS NEPHRO MALECOT è un catetere Malecot per drenaggio nefrostomico.

Caratteristiche del dispositivo:

- catetere in poliuretano, flessibile, radiopaco, con rivestimento idrofilico, resistente al kinking e all'ostruzione
- punta rinforzata per evitare l'incamiciamento durante l'inserimento
- punta Malecot
- connettore luer lock
- include mandrino e cannula di irrigidimento metallica a punta smussa per il modello per puntura diretta, cannula flessibile e cannula di irrigidimento metallica a punta smussa per il modello per tecnica Seldinger
- eccellente capacità di drenaggio.

HS NEPHRO MALECOT a nephrostomy drainage catheter.

Device Features:

- polyurethane catheter, flexible, radiopaque, with hydrophilic coating, kink and encrustation resistant
- reinforced tip to help prevent catheter from buckling during insertion
- Malecot type tip
- luer lock connector
- includes trocar and metal cannula (direct puncture) or flexible and metal cannula (Seldinger technique)
- excellent flow rates.

CODE	DIAMETER	WORKING LENGTH	BOX
Direct puncture			
HSNM0630	Fr. 6	290mm	5
HSNM0730	Fr. 7	290mm	5
HSNM0830	Fr. 8	290mm	5
HSNM1030	Fr. 10	290mm	5
HSNM1230	Fr. 12	290mm	5
HSNM1430	Fr. 14	290mm	5
HSNM1630	Fr. 16	290mm	5

CODE	DIAMETER	WORKING LENGTH	BOX
Seldinger technique			
HSNM0630S	Fr. 6	290mm	5
HSNM0730S	Fr. 7	290mm	5
HSNM0830S	Fr. 8	290mm	5
HSNM1030S	Fr. 10	290mm	5
HSNM1230S	Fr. 12	290mm	5
HSNM1430S	Fr. 14	290mm	5
HSNM1630S	Fr. 16	290mm	5

FUTUREMADEPRESENT

