

HS NOTA

DISPOSITIVO AUTOMATICO PER BIOPSIA ISTOLOGICA TRU-CUT AUTOMATIC DEVICE FOR TRU-CUT HISTOLOGICAL BIOPSY

HS NOTA è un ago a ghigliottina per biopsia dei tessuti molli.

HS NOTA is a soft tissue tru-cut biopsy needle.

Caratteristiche del dispositivo:

- leggero, maneggevole
- caricamento a rotazione
- indicatore del caricamento
- campionamento automatico
- dispositivo di sicurezza attivo con attivazione automatica durante la carica
- pulsante di sparo posteriore
- avanzamento 25mm o 15 mm
- marker ecogeno interno
- ago introduttore coassiale opzionale incluso nella stessa confezione (aggiungere "IN" alla fine del codice dell'ago)

Device Features:

- light, easy handling
- loading with rotating system
- loading indicator
- automatic sampling
- security feature that automatically activates during the device charging
- rear firing position
- 25mm or 15 mm penetration depth
- inner echogenic marker
- optional coaxial introducer needle included in the same primary packaging (add "IN" at the end of the needle code)

Avanzamento 25mm/ 25mm penetration depth


CODE	DIAMETER	LENGTH	BOX
HS NOTA14100	14G	100mm	10
HS NOTA14150	14G	150mm	10
HS NOTA14200	14G	200mm	10
HS NOTA16100	16G	100mm	10
HS NOTA16150	16G	150mm	10
HS NOTA16200	16G	200mm	10
HS NOTA16250	16G	250mm	10
HS NOTA18100	18G	100mm	10
HS NOTA18150	18G	150mm	10
HS NOTA18200	18G	200mm	10
HS NOTA18250	18G	250mm	10

Avanzamento 15mm/ 15mm penetration depth

CODE	DIAMETER	LENGTH	BOX
HS NOTA14100-1	14G	100mm	10
HS NOTA14150-1	14G	150mm	10
HS NOTA14200-1	14G	200mm	10
HS NOTA16100-1	16G	100mm	10
HS NOTA16150-1	16G	150mm	10
HS NOTA16200-1	16G	200mm	10
HS NOTA16250-1	16G	250mm	10
HS NOTA18100-1	18G	100mm	10
HS NOTA18150-1	18G	150mm	10
HS NOTA18200-1	18G	200mm	10
HS NOTA18250-1	18G	250mm	10



LEVA DI CARICAMENTO
TRIGGERING LEVER



FERMO DI SICUREZZA
SAFETY BLOCK



TRU-CUT Type

