

DISPOSITIVO SEMIAUTOMATICO TRU-CUT PER BIOPSIA ISTOLOGICA

TRU-CUT SEMIAUTOMATIC DEVICE FOR HISTOLOGICAL BIOPSY

PRECISA è un ago a ghigliottina per biopsia dei tessuti molli.

PRECISA is a soft tissue tru-cut biopsy needle.

Caratteristiche del dispositivo:


- cannula centimetrata sganciabile con cono luer lock, per utilizzo sotto TC e prelievi multipli
- mandrino rimovibile con intaglio da 20mm
- mandrino aggiuntivo a punta piramidale per il posizionamento sotto TC (non disponibile per la misura 13G)
- impugnatura ad anelli con finestra per la visualizzazione della dimensione del campione (10 o 20mm)
- marker ecogeno interno
- supporto per stereotassi
- fermo di sicurezza

Device Features:


- *releasable cannula with depth marks and luer lock connector for multiple biopsies and biopsy under CT guidance*
- *20mm sample notch removable stylet*
- *additional pyramidal tip stylet for positioning under CT guidance (not available for 13G)*
- *two rings grip; visualization of selected sample dimension (10 or 20 mm) through a dedicated aperture*
- *inner echogenic marker*
- *stereotactic spacer*
- *safety block*

CODE	DIAMETER	LENGTH	BOX
PRE1307	13G	70mm	10
PRE1310	13G	100mm	10
PRE1407	14G	70mm	10
PRE1410	14G	100mm	10
PRE1415	14G	150mm	10
PRE1420	14G	200mm	10
PRE1610	16G	100mm	10
PRE1615	16G	150mm	10
PRE1620	16G	200mm	10

CODE	DIAMETER	LENGTH	BOX
PRE1807	18G	70mm	10
PRE1810	18G	100mm	10
PRE1815	18G	150mm	10
PRE1820	18G	200mm	10
PRE1825	18G	250mm	10
PRE2010	20G	100mm	10
PRE2015	20G	150mm	10
PRE2020	20G	200mm	10



FERMO DI SICUREZZA BREVETTATO
PATENTED SAFETY BLOCK



CANNULA ESTERNA RIMOVIBILE
REMOVABLE CANNULA



TRU-CUT Type

