

# TRUE-CUT

## DISPOSITIVO SEMIAUTOMATICO TRU-CUT PER BIOPSIA ISTOLOGICA TRU-CUT SEMIAUTOMATIC DEVICE FOR HISTOLOGICAL BIOPSY

**TRUE-CUT** è un ago a ghigliottina per biopsia dei tessuti molli.

### Caratteristiche del dispositivo:

- cannula centimetrata con fermo di profondità
- mandrino con intaglio da 18mm
- possibilità di scegliere la dimensione del campione
- impugnatura ad anelli con finestra per la visualizzazione della dimensione del campione (10 o 20mm)
- marker ecogeno interno
- fermo di sicurezza
- ago introduttore coassiale opzionale (aggiungere "-I" alla fine del codice dell'ago)

**TRUE-CUT** is a soft tissue tru-cut biopsy needle.

### Device Features:

- cannula with depth marks and depth stop
- stylet with 18mm sample notch
- sample dimension choice possibility
- two rings grip with aperture to see the chosen sample dimension (10 o 20mm)
- inner echogenic marker
- safety block
- optional coaxial introducer needle (add "-I" at the end of the needle code)

| CODE   | DIAMETER | LENGTH | BOX |
|--------|----------|--------|-----|
| TC1410 | 14G      | 100mm  | 10  |
| TC1415 | 14G      | 150mm  | 10  |
| TC1420 | 14G      | 200mm  | 10  |
| TC1610 | 16G      | 100mm  | 10  |
| TC1615 | 16G      | 150mm  | 10  |
| TC1620 | 16G      | 200mm  | 10  |

| CODE   | DIAMETER | LENGTH | BOX |
|--------|----------|--------|-----|
| TC1810 | 18G      | 100mm  | 10  |
| TC1815 | 18G      | 150mm  | 10  |
| TC1820 | 18G      | 200mm  | 10  |
| TC2010 | 20G      | 100mm  | 10  |
| TC2015 | 20G      | 150mm  | 10  |
| TC2020 | 20G      | 200mm  | 10  |



TRU-CUT Type

