

BIOMOL KIT

KIT CON DISPOSITIVO SEMIAUTOMATICO
PER BIOPSIA ISTOLOGICA CON TECNICA MENGHINI

*KIT WITH SEMIAUTOMATIC DEVICE
FOR MENGHINI-TYPE BIOPSY*



BIOMOL KIT

KIT CON DISPOSITIVO SEMIAUTOMATICO PER BIOPSIA ISTOLOGICA CON TECNICA MENGHINI
KIT WITH SEMIAUTOMATIC DEVICE FOR MENGHINI-TYPE BIOPSY

BIOMOL KIT è un kit completo di ago per biopsia dei tessuti molli che automatizza la procedura, lasciando allo stesso tempo all'operatore la sensibilità del prelievo. Il BIOMOL permette di eseguire la procedura bioptica con una sola mano, rendendola più veloce.

Caratteristiche del dispositivo:

- cannula sganciabile (rende possibile la biopsia sotto TC ed il doppio prelievo) con punta Menghini, centimetrica, con fermo di profondità, a pareti ultrasottili.
- mandrino a punta piramidale solidale al pistone della siringa.
- siringa con blocco della depressione.
- marker ecogeno interno.
- evita sanguinamento e seeding

Il BIOMOL KIT contiene i seguenti accessori:

- composizione BM####K-00:
telo per campo sterile 60x60cm (2), provetta con formalina 10% 9ml, garze 7.5x7.5cm (10), cerotto 10x8cm, mini bisturi, applicatore con spugna, siringa da 10ml con aghi ipodermici da 18Gx40mm (pre-filtrato) e 21Gx40mm, lidocaina 2% 10ml.
- composizione BM####K-01:
telo per campo sterile 60x60cm (2), provetta con formalina 10% 9ml, garze 7.5x7.5cm (10), cerotto 10x8cm, mini bisturi, applicatore con spugna.

BIOMOL KIT is a soft tissues biopsy needle kit. It automates the biopsy procedure and at the same time it leaves to the physician full control the sampling.

The BIOMOL allows for faster biopsies with one hand only.

Needle features:

- removable extra-thin-walled needle (it allows biopsies under CT guidance and double sampling) with Menghini tip, depth marks and depth stop.
- plunger integrated stylet with pyramidal tip.
- depression block.
- inner echogenic marker.

The BIOMOL KIT includes:

- BM####K-00 accessories:
drape 60x60cm (2), formalin 10% 9ml, gauzes 7.5x7.5cm (10), adhesive bandage 10x8cm, scalpel, sponge applicator, 10ml syringe with hypodermic needle 18Gx40mm (pre-filtered) and 21Gx40mm, lidocaine 2% 10ml.
- BM####K-01 accessories:
drape 60x60cm (2), formalin 10% 9ml, gauzes 7.5x7.5cm (10), adhesive bandage 10x8cm, scalpel, sponge applicator.

CODEDIAMETER	LENGTH	BOX	
BM1607	16G	70mm	5
BM1609	16G	90mm	5
BM1612	16G	120mm	5
BM1615	16G	150mm	5
BM1620	16G	200mm	5
BM1630	16G	300mm	5
BM1709	17G	90mm	5
BM1715	17G	150mm	5
BM1720	17G	200mm	5
BM1807	17G	70mm	5
BM1810	18G	100mm	5
BM1815	18G	150mm	5
BM1820	18G	200mm	5
BM1910	19G	100mm	5

CODEDIAMETER	LENGTH	BOX	
BM1910	19G	100mm	5
BM1915	19G	150mm	5
BM1920	19G	200mm	5
BM2010	20G	100mm	5
BM2015	20G	150mm	5
BM2020	20G	200mm	5
BM2110	21G	100mm	5
BM2115	21G	150mm	5
BM2120	21G	200mm	5
BM2123	21G	230mm	5
BM2205	22G	50mm	5
BM2207	22G	70mm	5
BM2210	22G	100mm	5
BM2215	22G	150mm	5
BM2220	22G	200mm	5

Codice Kit: aggiungere K-00 o K-01 alla fine del codice dell'ago
Kit code: add K-00 or K-01 at the end of needle code