

ECOJEKT

AGO PER BIOPSIA CITOLOGICA
ED ALTRE APPLICAZIONI

*NEEDLE FOR FNAC AND
OTHER APPLICATIONS*

FUTUREMADEPRESENT



ECOJEKT

AGO PER BIOPSIA CITOLOGICA ED ALTRE APPLICAZIONI
NEEDLE FOR FNAC AND OTHER APPLICATIONS

ECOJEKT è un ago universale con punta Chiba.

Caratteristiche del dispositivo:

- ago ultra sottile con punta a tripla affilatura, centimetrato, con fermo di profondità.
- mandrino rimovibile con codice colore.
- marker ecogeno interno (non presente sulle misure 25G).
- applicazioni: biopsia citologica, agoaspirazione, amniocentesi, villocentesi, prelievo di ovociti, prelievo di sangue fetale, alcolizzazione, puntura esplorativa, iniezione di mezzi di contrasto o altre sostanze.

ECOJEKT is a universal needle with Chiba tip.

Device Features:

- extra-thin-walled needle with triple sharpening tip, depth marks and depth stop.
- removable colour coded stylet.
- inner echogenic marker (absent on 25G items).
- intended uses: cytological biopsy, fine-needle aspiration, amniocentesis, chorionic villus sampling, oocyte pick-up, fetal blood sampling, percutaneous, alcoholization, exploratory puncture, contrast medium and other fluids injection.

| CODE | DIAMETER | LENGTH | BOX |
|-------|----------|--------|-----|
| E1715 | 17G | 150mm | 10 |
| E1720 | 17G | 200mm | 10 |
| E1815 | 18G | 150mm | 10 |
| E1818 | 18G | 180mm | 10 |
| E1820 | 18G | 200mm | 10 |
| E1825 | 18G | 250mm | 10 |
| E1828 | 18G | 280mm | 10 |
| E1907 | 19G | 70mm | 10 |
| E1910 | 19G | 100mm | 10 |
| E1915 | 19G | 150mm | 10 |
| E1918 | 19G | 180mm | 10 |
| E1920 | 19G | 200mm | 10 |
| E1930 | 19G | 300mm | 10 |
| E2005 | 20G | 50mm | 10 |
| E2007 | 20G | 70mm | 10 |
| E2010 | 20G | 100mm | 10 |
| E2015 | 20G | 150mm | 10 |
| E2018 | 20G | 180mm | 10 |
| E2020 | 20G | 200mm | 10 |
| E2030 | 20G | 300mm | 10 |
| E2035 | 20G | 350mm | 10 |

| CODE | DIAMETER | LENGTH | BOX |
|--------|----------|--------|-----|
| E2107 | 21G | 70mm | 10 |
| E2110 | 21G | 100mm | 10 |
| E2115 | 21G | 150mm | 10 |
| E2118 | 21G | 180mm | 10 |
| E2120 | 21G | 200mm | 10 |
| E2130 | 21G | 300mm | 10 |
| E2205 | 22G | 50mm | 10 |
| E22075 | 22G | 75mm | 10 |
| E2210 | 22G | 100mm | 10 |
| E2215 | 22G | 150mm | 10 |
| E2218 | 22G | 180mm | 10 |
| E2220 | 22G | 200mm | 10 |
| E2315 | 23G | 150mm | 10 |
| E2318 | 23G | 180mm | 10 |
| E2320 | 23G | 200mm | 10 |
| E2507 | 25G | 70mm | 10 |
| E2510 | 25G | 100mm | 10 |
| E2515 | 25G | 150mm | 10 |