

HARVEST NEEDLE

AGO PER ESPIANTO DEL MIDOLLO OSSEO

BONE MARROW TRANSPLANT NEEDLE



FORI LATERALI
LATERAL HOLES



HARVEST NEEDLE

AGO PER ESPIANTO DEL MIDOLLO OSSEO
BONE MARROW TRANSPLANT NEEDLE

HARVEST NEEDLE è un ago per espianto con fori laterali.

Caratteristiche del dispositivo:

- ago con punta a corona per una facile penetrazione del tessuto osseo, con impugnatura ergonomica che consente di distribuire la forza su tutto il palmo della mano e non localmente, con cinque ampi fori laterali in prossimità della punta.
- connettore luer lock per una connessione sicura della siringa.
- mandrino estraibile, perfettamente accoppiato all'ago, con punta piramidale.
- mandrino di ricambio.

HARVEST NEEDLE is a bone marrow transplant needle with lateral holes.

Device Features:

- *needle with crown tip to facilitate bone penetration, with ergonomic grip (it ensures safe grip and manoeuvrability and force distribution on the palm of the hand) and five holes.*
- *luer lock connector on the handle for secure syringe attachment.*
- *removable stylet with pyramidal tip.*
- *replacement stylet*

CODE	DIAMETER	LENGTH	BOX
ES 1110.PD	11G	100mm	10
ES 1310.PD	13G	100mm	10