



HS GENERAL PURPOSE

CATETERE PER DRENAGGIO PERCUTANEO

PERCUTANEOUS DRAINAGE CATHETER



HS GENERAL PURPOSE

CATETERE PER DRENAGGIO PERCUTANEO

PERCUTANEOUS DRAINAGE CATHETER

HS GENERAL PURPOSE è un catetere per uso generale con tecnica diretta o Seldinger.

Caratteristiche del dispositivo:

- catetere realizzato in poliuretano, flessibile, radiopaco, con rivestimento idrofilico, resistente al kinking e all'ostruzione
- particolare disegno del lume interno che evita l'incamiciamento durante l'inserimento e facilita il ritiro della cannula di irrigidimento
- punta pigtail (2cm) rastremata con ampi fori
- "locking system" per il blocco del pigtail
- connettore luer lock
- raddrizzatore per la punta, mandrino a punta piramidale e cannula di irrigidimento metallica a punta smussa
- eccellente capacità di drenaggio.

HS GENERAL PURPOSE is a general purpose catheter intended to be used with one-shot or Seldinger technique.

Device Features:

- polyurethane catheter, flexible, radiopaque, with hydrophilic coating, kink and encrustation resistant
- reinforced tip to help prevent catheter from buckling during insertion
- 2cm pigtail tip with large holes
- locking mechanism for pigtail tip
- luer lock connector
- includes split straightener, trocar and metal cannula
- excellent flow rates.

| CODE | DIAMETER | LENGTH | BOX |
|---|----------|--------|-----|
| Cateteri senza locking system/Non-locking catheters | | | |
| HSGP0625 | Fr. 6 | 240mm | 5 |
| HSGP0630 | Fr. 6 | 290mm | 5 |
| HSGP0725 | Fr. 7 | 240mm | 5 |
| HSGP0730 | Fr. 7 | 290mm | 5 |
| HSGP0740 | Fr. 7 | 390mm | 5 |
| HSGP0825 | Fr. 8 | 240mm | 5 |
| HSGP0830 | Fr. 8 | 290mm | 5 |
| HSGP0840 | Fr. 8 | 390mm | 5 |
| HSGP1030 | Fr. 10 | 290mm | 5 |
| HSGP1040 | Fr. 10 | 390mm | 5 |
| HSGP1230 | Fr. 12 | 290mm | 5 |
| HSGP1240 | Fr. 12 | 390mm | 5 |
| HSGP1430 | Fr. 14 | 290mm | 5 |
| HSGP1440 | Fr. 14 | 390mm | 5 |
| HSGP1630 | Fr. 16 | 290mm | 5 |
| HSGP1640 | Fr. 16 | 390mm | 5 |

| CODE | DIAMETER | LENGTH | BOX |
|--|----------|--------|-----|
| Cateteri con locking system /Locking catheters | | | |
| HSGPL0625 | Fr. 6 | 240mm | 5 |
| HSGPL0630 | Fr. 6 | 290mm | 5 |
| HSGPL0725 | Fr. 7 | 240mm | 5 |
| HSGPL0730 | Fr. 7 | 290mm | 5 |
| HSGPL0740 | Fr. 7 | 390mm | 5 |
| HSGPL0825 | Fr. 8 | 240mm | 5 |
| HSGPL0830 | Fr. 8 | 290mm | 5 |
| HSGPL0840 | Fr. 8 | 390mm | 5 |
| HSGPL1030 | Fr. 10 | 290mm | 5 |
| HSGPL1040 | Fr. 10 | 390mm | 5 |
| HSGPL1230 | Fr. 12 | 290mm | 5 |
| HSGPL1240 | Fr. 12 | 390mm | 5 |
| HSGPL1430 | Fr. 14 | 290mm | 5 |
| HSGPL1440 | Fr. 14 | 390mm | 5 |
| HSGPL1630 | Fr. 16 | 290mm | 5 |
| HSGPL1640 | Fr. 16 | 390mm | 5 |

* lunghezza operativa/working length

I cateteri Fr.6-7 accettano guide 0.035"/Fr.6-7 catheters accept guidewire 0.035".
I cateteri Fr.8-16 accettano guide 0.038"/Fr.8-16 catheters accept guidewire 0.038".