

# HS GENERAL PURPOSE KIT

KIT CON CATETERE PER DRENAGGIO PERCUTANEO

*PERCUTANEOUS DRAINAGE CATHETER KIT*

**HS GENERAL PURPOSE KIT** è un catetere per uso generale con tecnica diretta o Seldinger.

## Caratteristiche del dispositivo:

- catetere realizzato in poliuretano, flessibile, radiopaco, con rivestimento idrofilico, resistente al kinking e all'ostruzione
- particolare disegno del lume interno che evita l'incamiciamento durante l'inserimento e facilita il ritiro della cannula di irrigidimento
- punta pigtail (2cm) rastremata con ampi fori
- "locking system" per il blocco del pigtail
- connettore luer lock
- raddrizzatore per la punta, mandrino a punta piramidale e cannula di irrigidimento metallica a punta smussa
- eccellente capacità di drenaggio.

**HS GENERAL PURPOSE KIT** contiene i seguenti accessori: ago introduttore 18G, guida flessibile o tipo Lunderquist (teflonata, idrofilica, con tratto terminale floppy e punta a "J" 3mm), dilatatori sequenziali con punta rastremata, raccordo di estensione con connettori luer lock, raccordo per sacca con connettori luer lock e ad imbuto, ago Ecojekt 22G, disco di stabilizzazione in silicone, valvola per il lavaggio del catetere, telo per campo sterile 60x60cm, salvietta bianca 40x40cm, garze 7.5x7.5cm (5), rubinetto a tre vie, bisturi, applicatore con spugna, pinzetta, gel sterile.

**HS GENERAL PURPOSE KIT** is a general purpose catheter intended to be used with one-shot or Seldinger technique.

## Device Features:

- tri-axial system with inner trapping device
- throw length 25mm or 15mm
- light, easy handling, easy to use and to charge
- loading with retraction system
- activation indicator
- automatic sampling
- lateral or rear firing button
- passive safety system
- safe for patients
- avoids bleeding and seeding

## **HS GENERAL PURPOSE KIT** includes:

introducer needle 18G, flexible or Lunderquist type guidewire (PTFE-coated, hydrophilic, with floppy and 3mm "J" tip), dilators, luer lock extension tube, funnel-luer lok adapter, Ecojekt needle 22G, catheter stabilizer, septum cap, drape 60x60cm, white drape 40x40cm, gauzes 7.5x7.5cm (5), three-way stopcock, scalpel, sponge applicator, forceps, sterile gel.

# HS GENERAL PURPOSE KIT

KIT CON CATETERE PER DRENAGGIO PERCUTANEO  
PERCUTANEOUS DRAINAGE CATHETER KIT

CODE	DIAMETER	LENGTH	BOX
<i>Cateteri senza locking system e guida flessibile/ Non-locking catheters and flexible guidewire</i>			
HSGP0625KF	Fr. 6	240mm	1
HSGP0630KF	Fr. 6	290mm	1
HSGP0725KF	Fr. 7	240mm	1
HSGP0730KF	Fr. 7	290mm	1
HSGP0740KF	Fr. 7	390mm	1
HSGP0825KF	Fr. 8	240mm	1
HSGP0830KF	Fr. 8	290mm	1
HSGP0840KF	Fr. 8	390mm	1
HSGP1030KF	Fr. 10	290mm	1
HSGP1040KF	Fr. 10	390mm	1
HSGP1230KF	Fr. 12	290mm	1
HSGP1240KF	Fr. 12	390mm	1
HSGP1430KF	Fr. 14	290mm	1
HSGP1440KF	Fr. 14	390mm	1
HSGP1630KF	Fr. 16	290mm	1
HSGP1640KF	Fr. 16	390mm	1

CODE	DIAMETER	LENGTH	BOX
<i>Cateteri con locking system e guida flessibile/ Locking catheters and flexible guidewire</i>			
HSGPL0625KF	Fr. 6	240mm	1
HSGPL0630KF	Fr. 6	290mm	1
HSGPL0725KF	Fr. 7	240mm	1
HSGPL0730KF	Fr. 7	290mm	1
HSGPL0740KF	Fr. 7	390mm	1
HSGPL0825KF	Fr. 8	240mm	1
HSGPL0830KF	Fr. 8	290mm	1
HSGPL0840KF	Fr. 8	390mm	1
HSGPL1030KF	Fr. 10	290mm	1
HSGPL1040KF	Fr. 10	390mm	1
HSGPL1230KF	Fr. 12	290mm	1
HSGPL1240KF	Fr. 12	390mm	1
HSGPL1430KF	Fr. 14	290mm	1
HSGPL1440KF	Fr. 14	390mm	1
HSGPL1630KF	Fr. 16	290mm	1
HSGPL1640KF	Fr. 16	390mm	1

CODE	DIAMETER	LENGTH	BOX
<i>Cateteri senza locking system e guida tipo Lunderquist/ Non-locking catheters and Lunderquist guidewire</i>			
HSGP0625KL	Fr. 6	240mm	1
HSGP0630KL	Fr. 6	290mm	1
HSGP0725KL	Fr. 7	240mm	1
HSGP0730KL	Fr. 7	290mm	1
HSGP0740KL	Fr. 7	390mm	1
HSGP0825KL	Fr. 8	240mm	1
HSGP0830KL	Fr. 8	290mm	1
HSGP0840KL	Fr. 8	390mm	1
HSGP1030KL	Fr. 10	290mm	1
HSGP1040KL	Fr. 10	390mm	1
HSGP1230KL	Fr. 12	290mm	1
HSGP1240KL	Fr. 12	390mm	1
HSGP1430KL	Fr. 14	290mm	1
HSGP1440KL	Fr. 14	390mm	1
HSGP1630KL	Fr. 16	290mm	1
HSGP1640KL	Fr. 16	390mm	1

CODE	DIAMETER	LENGTH	BOX
<i>Cateteri con locking system e guida tipo Lunderquist/ Locking catheters and Lunderquist guidewire</i>			
HSGPL0625KL	Fr. 6	240mm	1
HSGPL0630KL	Fr. 6	290mm	1
HSGPL0725KL	Fr. 7	240mm	1
HSGPL0730KL	Fr. 7	290mm	1
HSGPL0740KL	Fr. 7	390mm	1
HSGPL0825KL	Fr. 8	240mm	1
HSGPL0830KL	Fr. 8	290mm	1
HSGPL0840KL	Fr. 8	390mm	1
HSGPL1030KL	Fr. 10	290mm	1
HSGPL1040KL	Fr. 10	390mm	1
HSGPL1230KL	Fr. 12	290mm	1
HSGPL1240KL	Fr. 12	390mm	1
HSGPL1430KL	Fr. 14	290mm	1
HSGPL1440KL	Fr. 14	390mm	1
HSGPL1630KL	Fr. 16	290mm	1
HSGPL1640KL	Fr. 16	390mm	1

\* lunghezza operativa/working length