

HS GENERAL PURPOSE STRAIGHT

CATETERE PER DRENAGGIO PERCUTANEO

PERCUTANEOUS DRAINAGE CATHETER

FUTUREMADEPRESENT



HS GENERAL PURPOSE STRAIGHT

CATETERE PER DRENAGGIO PERCUTANEO
PERCUTANEOUS DRAINAGE CATHETER

HS GENERAL PURPOSE STRAIGHT è un catetere per uso generale con tecnica diretta.

Caratteristiche del dispositivo:

- catetere realizzato in poliuretano, flessibile, radiopaco, con rivestimento idrofilico, resistente al kinking e all'ostruzione
- particolare disegno del lume interno che evita l'incamiciamento durante l'inserimento e facilita il ritiro della cannula di irrigidimento
- punta retta rastremata con ampi fori
- connettore luer lock
- mandrino a punta piramidale e cannula di irrigidimento metallica a punta smussa
- eccellente capacità di drenaggio

HS GENERAL PURPOSE STRAIGHT is a general purpose catheter intended to be used with one-shot technique.

Device Features:

- polyurethane catheter, flexible, radiopaque, with hydrophilic coating, kink and encrustation resistant
- reinforced tip to help prevent catheter from buckling during insertion
- straight tip with large holes
- luer lock connector
- includes trocar and metal cannula
- excellent flow rates.

CODE	DIAMETER	LENGTH	BOX
HSGPS0620	Fr. 6	190mm	5
HSGPS0825	Fr. 8	240mm	5

* lunghezza operativa/working length