

HS NEPHRO CATHETER KIT

KIT CON CATETERE PER DRENAGGIO NEFROSTOMICO

NEPHROSTOMY DRAINAGE CATHETER KIT

HS NEPHRO CATHETER KIT è un kit con catetere per drenaggio nefrostomico.

Caratteristiche del dispositivo:

- catetere realizzato in poliuretano, flessibile, radiopaco, con rivestimento idrofilico, resistente al kinking e all'ostruzione
- particolare disegno del lume interno che evita l'incamiciamento durante l'inserimento e facilita il ritiro della cannula di irrigidimento
- punta pigtail (2cm) rastremata con ampi fori
- "locking system" per il blocco del pigtail
- connettore luer lock
- raddrizzatore per la punta, cannula di irrigidimento metallica a punta smussa e cannula flessibile
- eccellente capacità di drenaggio.

HS NEPHRO CATHETER KIT contiene i seguenti accessori: ago introduttore 18G, guida flessibile o tipo Lunderquist (teflonata, idrofilica, con tratto terminale floppy e punta a "J" 3mm), dilatatori sequenziali con punta rastremata, raccordo di estensione con connettori luer lock, raccordo per sacca con connettori luer lock e ad imbuto, ago Ecojekt 22G, disco di stabilizzazione in silicone, valvola per il lavaggio del catetere, telo per campo sterile 60x60cm, salvietta bianca 40x40cm, garze 7.5x7.5cm (5), rubinetto a tre vie, bisturi, applicatore con spugna, pinzetta, gel sterile.

HS NEPHRO CATHETER KIT is a nephrostomy drainage catheter.

Device Features:

- polyurethane catheter, flexible, radiopaque, with hydrophilic coating, kink and encrustation resistant
- reinforced tip to help prevent catheter from buckling during insertion
- 2cm pigtail tip with large holes
- locking mechanism for pigtail tip
- luer lock connector
- includes split straightener, metal cannula and flexible cannula
- excellent flow rates.

HS NEPHRO CATHETER KIT includes:

introducer needle 18G, flexible or Lunderquist type guidewire (PTFE-coated, hydrophilic, with floppy and 3mm "J" tip), dilators, luer lock extension tube, funnel-luer lok adapter, Ecojekt needle 22G, catheter stabilizer, septum cap, drape 60x60cm, white drape 40x40cm, gauzes 7.5x7.5cm (5), three-way stopcock, scalpel, sponge applicator, forceps, sterile gel.

HS NEPHRO CATHETER KIT

KIT CON CATETERE PER DRENAGGIO NEFROSTOMICO
NEPHROSTOMY DRAINAGE CATHETER KIT

CODE	DIAMETER	LENGTH	BOX
<i>Cateteri senza locking system e guida flessibile/ Non-locking catheters and flexible guidewire</i>			
HSNP0740KF	Fr. 7	390mm	1
HSNP0750KF	Fr. 7	490mm	1
HSNP0840KF	Fr. 8	390mm	1
HSNP0850KF	Fr. 8	490mm	1
HSNP1040KF	Fr. 10	390mm	1
HSNP1050KF	Fr. 10	490mm	1
HSNP1240KF	Fr. 12	390mm	1
HSNP1250KF	Fr. 12	490mm	1
HSNP1440KF	Fr. 14	390mm	1
HSNP1450KF	Fr. 14	490mm	1

CODE	DIAMETER	LENGTH	BOX
<i>Cateteri con locking system e guida flessibile/ Locking catheters and flexible guidewire</i>			
HSNPL0740KF	Fr. 7	390mm	1
HSNPL0750KF	Fr. 7	490mm	1
HSNPL0840KF	Fr. 8	390mm	1
HSNPL0850KF	Fr. 8	490mm	1
HSNPL1040KF	Fr. 10	390mm	1
HSNPL1050KF	Fr. 10	490mm	1
HSNPL1240KF	Fr. 12	390mm	1
HSNPL1250KF	Fr. 12	490mm	1
HSNPL1440KF	Fr. 14	390mm	1
HSNPL1450KF	Fr. 14	490mm	1

CODE	DIAMETER	LENGTH	BOX
<i>Cateteri senza locking system e guida tipo Lunderquist/ Non-locking catheters and Lunderquist guidewire</i>			
HSNP0740KL	Fr. 7	390mm	1
HSNP0750KL	Fr. 7	490mm	1
HSNP0840KL	Fr. 8	390mm	1
HSNP0850KL	Fr. 8	490mm	1
HSNP1040KL	Fr. 10	390mm	1
HSNP1050KL	Fr. 10	490mm	1
HSNP1240KL	Fr. 12	390mm	1
HSNP1250KL	Fr. 12	490mm	1
HSNP1440KL	Fr. 14	390mm	1
HSNP1450KL	Fr. 14	490mm	1

CODE	DIAMETER	LENGTH	BOX
<i>Cateteri con locking system e guida tipo Lunderquist/ Locking catheters and Lunderquist guidewire</i>			
HSNPL0740KL	Fr. 7	390mm	1
HSNPL0750KL	Fr. 7	490mm	1
HSNPL0840KL	Fr. 8	390mm	1
HSNPL0850KL	Fr. 8	490mm	1
HSNPL1040KL	Fr. 10	390mm	1
HSNPL1050KL	Fr. 10	490mm	1
HSNPL1240KL	Fr. 12	390mm	1
HSNPL1250KL	Fr. 12	490mm	1
HSNPL1440KL	Fr. 14	390mm	1
HSNPL1450KL	Fr. 14	490mm	1

* lunghezza operativa/working length

I cateteri Fr.6-7 accettano guide 0.035"/Fr.6-7 catheters accept guidewire 0.035".
I cateteri Fr.8-14 accettano guide 0.038"/Fr.8-14 catheters accept guidewire 0.038".