


INTRO-KIT

SET PER L'INTRODUZIONE
DI CATERI PERCUTANEI

INTRODUCER SET



INTRO-KIT

SET PER L'INTRODUZIONE DI CATETERI PERCUTANEI

INTRODUCER SET

INTRO-KIT è un set per l'introduzione dei cateteri percutanei con tecnica Seldinger.

Caratteristiche del kit:

- ago introduttore Ecojekt BS 18G.
- guida flessibile o tipo Lunderquist (teflonata, idrofilica, con tratto terminale floppy e punta a "J" 3mm).
- dilatatori sequenziali con punta rastremata.
- raccordo di estensione con connettori luer lock.
- raccordo per sacca con connettori luer lock e ad imbuto.
- ago Ecojekt 22G.
- disco di stabilizzazione in silicone.
- valvola per il lavaggio del catetere.

INTRO-KIT contiene i seguenti accessori:

telo per campo sterile 60x60cm, salvietta bianca 40x44cm, garze 7.5x7.5cm (5), rubinetto a tre vie, bisturi, applicatore con spugna, gel sterile.

INTRO-KIT provides for easy, accurate, non-traumatic over-the-guidewire catheter placement.

Kit features:

- introducer needle Ecojekt BS 18G.
- flexible or Lunderquist type guidewire (PTFE-coated, hydrophilic, with floppy and 3mm "J" tip).
- dilators.
- luer lock extension tube.
- funnel-luer lock adapter.
- Ecojekt needle 22G.
- catheter stabilizer.
- septum cap.

INTRO-KIT includes:

drape 60x60cm, white drape 40x44cm, gauzes 7.5x7.5cm (5), three-way stopcock, scalpel, sponge applicator, sterile gel.

CODE	DIAMETER	BOX
Kit con guida flessibile Kit with flexible guidewire 0.035"x1000mm		
INTRO8F	per catetere/for catheter Fr.8	5
INTRO10F	per catetere/for catheter Fr.10	5
INTRO12F	per catetere/for catheter Fr.12	5
INTRO14F	per catetere/for catheter Fr.14	5
Kit con guida tipo Lunderquist Kit with Lunderquist type guidewire 0.035"x1000mm		
INTRO8L	per catetere/for catheter Fr.8	5
INTRO10L	per catetere/for catheter Fr.10	5
INTRO12L	per catetere/for catheter Fr.12	5
INTRO14L	per catetere/for catheter Fr.14	5