

QUINJEKT

AGO PER BIOPSIA CITOLOGICA
ED ALTRE APPLICAZIONI

*NEEDLE FOR FNAC
AND OTHER APPLICATIONS*



QUINJEKT

AGO PER BIOPSIA CITOLOGICA ED ALTRE APPLICAZIONI
NEEDLE FOR FNAC AND OTHER APPLICATIONS

QUINJEKT® è un ago universale con punta tipo Quincke.

Caratteristiche del dispositivo:

- ago con punta spinale.
- mandrino rimovibile con codice colore.
- marker ecogeno interno.
- applicazioni: biopsia citologica, agoaspirazione, amniocentesi, villocentesi, prelievo di ovociti, prelievo di sangue fetale, alcolizzazione, puntura esplorativa, iniezione di mezzi di contrasto o altre sostanze.

QUINJEKT is a universal needle with Quincke type tip.

Device Features:

- spinal tip needle.
- removable colour coded stylet.
- inner echogenic marker.
- intended uses: cytological biopsy, fine-needle aspiration, amniocentesis, chorionic villus sampling, oocyte pick-up, fetal blood sampling, percutaneous, alcoholization, exploratory puncture, contrast medium and other fluids injection.

CODE	DIAMETER	LENGTH	BOX
Q1805	18G	50mm	10
Q1809	18G	90mm	10
Q1812	18G	120mm	10
Q1815	18G	150mm	10
Q1820	18G	200mm	10
Q1825	18G	250mm	10
Q2005	20G	50mm	10
Q2009	20G	90mm	10
Q2012	20G	120mm	10
Q2015	20G	150mm	10
Q2020	20G	200mm	10
Q2025	20G	250mm	10
Q2105	21G	50mm	10
Q2109	21G	90mm	10

CODE	DIAMETER	LENGTH	BOX
Q2115	21G	150mm	10
Q2207	22G	70mm	10
Q2209	22G	90mm	10
Q2212	22G	120mm	10
Q22125	22G	125mm	10
Q2215	22G	150mm	10
Q22175	22G	175mm	10
Q2220	22G	200mm	10
Q2305	23G	50mm	10
Q2309	23G	90mm	10
Q2315	23G	150mm	10
Q2320	23G	200mm	10
Q2505	25G	50mm	10
Q2515	25G	150mm	10