

# SUCTION RECTAL BIOPSY SYSTEM

DISPOSITIVO PER BIOPSIA RETTALE

SUCTION RECTAL BIOPSY DEVICE



VASCHETTA PER STERILIZZAZIONE  
STERILIZATION STORAGE CONTAINER



MANOMETRO  
SUCTION METER



CAPSULE  
CAPSULES

# SUCTION RECTAL BIOPSY SYSTEM

DISPOSITIVO PER BIOPSIA RETTALE

SUCTION RECTAL BIOPSY KIT

**SUCTION RECTAL BIOPSY SYSTEM** è indicato per la biopsia rettale e per la diagnosi del morbo di Hirschsprung

## Caratteristiche del dispositivo:

- disegnato appositamente per ottenere un campione di mucosa e sottomucosa in modo controllato e costante.
- pistola con impugnatura ergonomica, composta da 4 pezzi, facile da smontare e rimontare.
- capsule cilindriche monouso con lama incapsulata che consentono un recupero semplice e veloce del campione.
- raccordo di suzione.
- manometro.
- contenitore pluriuso per sterilizzare la pistola e i suoi componenti.

**SUCTION RECTAL BIOPSY SYSTEM** is intended for rectal biopsy and Hirschsprung's disease diagnosis.

## Device Features:

- for controlled, consistent rectal biopsies of mucosa and submucosa.
- 4 pieces reusable gun with ergonomic handle, easy assembly and disassembly.
- disposable cylindrical capsules with enclosed blade, for easy and fast sample retrieval.
- suction tube.
- suction meter.
- sterilize storage container.

CODE	DESCRIPTION	BOX
CP1200	SUCTION RECTAL BIOPSY KIT(3 capsules)	10
MR100	SUCTION METER	10
HP1000	SUCTION RECTAL BIOPSY GUN	1
SS1000	STORAGE STERILIZATION TRAY	1