


TOP CUT

AGO TRU-CUT PER BIOPSIA ISTOLOGICA

TRU-CUT NEEDLE FOR HISTOLOGICAL BIOPSY



TRU-CUT Type

TOP CUT

AGO TRU-CUT PER BIOPSIA ISTOLOGICA
TRU-CUT NEEDLE FOR HISTOLOGICAL BIOPSY

TOP CUT è un ago a ghigliottina per i dispositivi automatici riutilizzabili HS GUN 03, BB GUN e PRO-MAG.

Caratteristiche del dispositivo:

- ago centimetrato con codice colore.
- mandrino rimovibile con intaglio da 18mm.
- distanziatore per il corretto posizionamento dell'ago all'interno della pistola.
- avanzamento 22mm o 11mm.
- marker ecogeno interno.

TOP CUT is a tru-cut needle intended to be used with HS GUN03, BB GUN and PRO-MAG biopsy guns.

Device Features:

- colour coded needle with depth marks.
- 18mm sample notch removable stylet.
- spacer for easier positioning of the needle into the biopsy gun.
- 22mm or 11mm penetration depth.
- inner echogenic marker.

CODE	DIAMETER	LENGTH	BOX
TOP CUT (AVANZAMENTO 22mm) (PENETRATION DEPTH 22mm)			
TOP 1410	14G	100mm	10
TOP 1415	14G	150mm	10
TOP 1420	14G	200mm	10
TOP 1610	16G	100mm	10
TOP 1615	16G	150mm	10
TOP 1620	16G	200mm	10
TOP 1625	16G	250mm	10
TOP 1810	18G	100mm	10
TOP 1815	18G	150mm	10
TOP 1820	18G	200mm	10
TOP 1825	18G	250mm	10
TOP 1830	18G	300mm	10
TOP 2010	20G	100mm	10
TOP 2015	20G	150mm	10
TOP 2020	20G	200mm	10

CODE	DIAMETER	LENGTH	BOX
TOP CUT SHORT (AVANZAMENTO 11mm) (PENETRATION DEPTH 11mm)			
TOP 14075 S	14G	75mm	10
TOP 14115 S	14G	115mm	10
TOP 16075 S	16G	75mm	10
TOP 16115 S	16G	115mm	10
TOP 18075 S	18G	75mm	10
TOP 18115 S	18G	115mm	10

PRO-MAG™ è un marchio della Argon Medical
PRO-MAG™ is a trademark of Argon Medical