


TRUE-CUT

DISPOSITIVO SEMIAUTOMATICO
TRU-CUT PER BIOPSIA ISTOLOGICA

*TRU-CUT SEMIAUTOMATIC DEVICE
FOR HISTOLOGICAL BIOPSY*



TRU-CUT Type

TRUE-CUT

DISPOSITIVO SEMIAUTOMATICO TRU-CUT PER BIOPSIA ISTOLOGICA

TRU-CUT SEMIAUTOMATIC DEVICE FOR HISTOLOGICAL BIOPSY

TRUE-CUT è un ago a ghiagliottina per biopsia dei tessuti molli.

TRUE-CUT is a soft tissues tru-cut biopsy needle.

Caratteristiche del dispositivo:

- cannula centimetrata con fermo di profondità
- mandrino con intaglio da 18mm
- impugnatura ad anelli per effettuare la biopsia con una sola mano, con finestra per la visualizzazione della dimensione del campione scelta
- marker ecogeno interno
- fermo di sicurezza.

Device Features:

- cannula with depth marks and depth stop
- stylet with 18mm sample notch
- two rings grip to perform the procedure with one hand only, with window to see the chosen sample dimension
- inner echogenic marker
- safety block.

CODE	DIAMETER	LENGTH	BOX
TC1410	14G	100mm	10
TC1415	14G	150mm	10
TC1420	14G	200mm	10
TC1610	16G	100mm	10
TC1615	16G	150mm	10
TC1620	16G	200mm	10

CODE	DIAMETER	LENGTH	BOX
TC1810	18G	100mm	10
TC1815	18G	150mm	10
TC1820	18G	200mm	10
TC2010	20G	100mm	10
TC2015	20G	150mm	10
TC2020	20G	200mm	10